

Datasheet: AHP420T

BATCH NUMBER 1608

Description:	RABBIT ANTI LAMININ
Specificity:	LAMININ
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/200 - 1/1000
Immunohistology - Paraffin	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Species Cross Reactivity

Reacts with: Rabbit, Chicken, Duck, Birds, Human, Dog

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Whole serum - liquid

Antiserum Preparation

Antisera to laminin were raised by repeated immunisation of rabbits with highly purified antigen.

Preservative Stabilisers

0.09% Sodium Azide

Immunogen

Purified laminin from Engelbreth-Holm-Swarm (EHS) murine sarcoma.

RRID

AB_2133776

Specificity	Rabbit anti Laminin antibody recognizes both A (~400 kDa) and B (~220 kDa) chains of laminin.
Histology Positive Control Tissue	Human kidney
References	<ol style="list-style-type: none"> 1. Gordon-Weeks, P.R. <i>et al.</i> (1989) Transient expression of laminin immunoreactivity in the developing rat hippocampus. J Neurocytol. 18 (4): 451-63. 2. Gordon-Weeks, P.R. <i>et al.</i> (1992) A study of the expression of laminin in the spinal cord of the frog during development and regeneration. Exp Physiol. 77 (5): 681-92. 3. Abd-Elmaksoud, A. <i>et al.</i> (2009) Comparative expression of laminin and smooth muscle actin in the testis and epididymis of poultry and rabbit. J Mol Histol. 40: 407-16. 4. Au, E. <i>et al.</i> (2007) SPARC from olfactory ensheathing cells stimulates Schwann cells to promote neurite outgrowth and enhances spinal cord repair. J Neurosci. 27: 7208-21. 5. Mueller, S.N. <i>et al.</i> (2007) Viral targeting of fibroblastic reticular cells contributes to immunosuppression and persistence during chronic infection. Proc Natl Acad Sci U S A. 104: 15430-5. 6. Falk, V. <i>et al.</i> (2002) Regulation of matrix metalloproteinases and effect of MMP-inhibition in heart transplant related reperfusion injury. Eur J Cardiothorac Surg. 22: 53-8. 7. Richardson, G.D. <i>et al.</i> (2009) Dynamic expression of Syndecan-1 during hair follicle morphogenesis. Gene Expr Patterns. 9: 454-60. 8. Miners, J.S. <i>et al.</i> (2008) Angiotensin-converting enzyme (ACE) levels and activity in Alzheimer's disease, and relationship of perivascular ACE-1 to cerebral amyloid angiopathy. Neuropathol Appl Neurobiol. 34: 181-93. 9. Melo, F.A. <i>et al.</i> (2009) Hepatic extracellular matrix alterations in dogs naturally infected with Leishmania (Leishmania) chagasi. Int J Exp Pathol. 90: 538-48. 10. Yeh, T.S. <i>et al.</i> (2014) Baculovirus-transduced, VEGF-expressing adipose-derived stem cell sheet for the treatment of myocardium infarction. Biomaterials. 35 (1): 174-84. 11. Foditsch, E.E. <i>et al.</i> (2016) Skeletal muscle proteins: a new approach to delimitate the time since death. Int J Legal Med. 130 (2): 433-40. 12. Kato, Y. <i>et al.</i> (2015) Targeting Antigen to Clec9A Primes Follicular Th Cell Memory Responses Capable of Robust Recall. J Immunol. 195 (3): 1006-14.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10081 available at: https://www.bio-rad-antibodies.com/SDS/AHP420T 10081

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M364210:200529'

Printed on 27 Apr 2023
