

Datasheet: AHP440

BATCH NUMBER 146282

Description:	RABBIT ANTI HUMAN CD253 (C-TERMINAL)
Specificity:	CD253 (C-TERMINAL)
Other names:	TRAIL
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/500

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to CD253 (TRAIL) were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	Peptide from C-terminus of human CD253 (TRAIL).
External Database Links	<p>UniProt: P50591 Related reagents</p> <p>Entrez Gene: 8743 TNFSF10 Related reagents</p>
Synonyms	APO2L, TRAIL
RRID	AB_322739
Specificity	Rabbit anti Human CD253 antibody recognizes the ~32 kDa apoptosis related protein designated TRAIL (TNF-related apoptosis-inducing ligand) CD253 also known as TRAIL. CD253 is a type II membrane protein expressed in a variety of tissues. Receptors for TRAIL have been identified and termed DR4 and DR5 (Griffith and Lynch 1998).
Further Reading	<ol style="list-style-type: none"> 1. Wiley, S.R. <i>et al.</i> (1995) Identification and characterization of a new member of the TNF family that induces apoptosis. Immunity. 3 (6): 673-82. 2. Pitti, R.M. <i>et al.</i> (1996) Induction of apoptosis by Apo-2 ligand, a new member of the tumor necrosis factor cytokine family. J Biol Chem. 271 (22): 12687-90.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP440 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M364221:200529'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)