

## Datasheet: AHP479

**BATCH NUMBER 172220**

<b>Description:</b>	RABBIT ANTI GFR ALPHA 2
<b>Specificity:</b>	GFR ALPHA 2
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/1000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

#### Species Cross Reactivity

Reacts with: Rat, Human, Mouse  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antisera to human GFR alpha-2 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by ion exchange chromatography.

**Buffer Solution** Phosphate buffered saline

**Preservative** 0.02% Sodium Azide

## Stabilisers

---

**Approx. Protein Concentrations** IgG concentration 1 mg/ml

---

**Immunogen** Synthetic peptide from C-terminus of human GFR alpha-2.

---

## External Database Links

**UniProt:**

[O00451](#)   [Related reagents](#)

**Entrez Gene:**

[2675](#)   GFRA2   [Related reagents](#)

---

**Synonyms** GDNFRB, RETL2, TRNR2

---

**RRID** AB\_2247568

---

**Specificity** **Rabbit anti GFR alpha-2 antibody** recognizes the cell surface receptor for neurturin, designated GFR alpha-2. Neurturin is a member of the GDNF family, and has an important role in the control of vertebrate neuron survival and differentiation. GRF alpha-2 is a ~52 kDa glycoprotein, and is a GPI-linked cell surface receptor.

---

## Further Reading

1. Kotzbauer, P.T. *et al.* (1996) Neurturin, a relative of glial-cell-line-derived neurotrophic factor. [Nature. 384 \(6608\): 467-70.](#)
  2. Klein, R.D. *et al.* (1997) A GPI-linked protein that interacts with Ret to form a candidate neurturin receptor. [Nature. 387 \(6634\): 717-21.](#)
  3. Buj-bello, A. *et al.* (1997) Neurturin responsiveness requires a GPI-linked receptor and the Ret receptor tyrosine kinase. [Nature. 387 \(6634\): 721-4.](#)
  4. GFR $\alpha$  Nomenclature Committee (1997) Nomenclature of GPI-linked receptors for the GDNF ligand family. [Neuron 19: 485.](#)
- 

## Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP479>

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)  
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)  
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)  
Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)  
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

### **Recommended Useful Reagents**

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M382377:210513'

**Printed on 29 Jan 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)