

Datasheet: AHP486

BATCH NUMBER 170130

Description:	RABBIT ANTI MOUSE CAD (aa205-222)
Specificity:	CAD (aa205-222)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/500 - 1/1000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species Mouse

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to murine CAD were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.02% Sodium Azide

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Synthetic peptide corresponding to amino acids 205-222 of mouse CAD.
External Database Links	<p>UniProt: O54788 Related reagents</p> <p>Entrez Gene: 13368 Dffb Related reagents</p>
Synonyms	Cad
RRID	AB_322892
Specificity	<p>Rabbit anti Murine CAD antibody recognizes murine CAD (Caspase Activated Deoxyribonuclease), the murine homologue of human DFF-40.</p> <p>CAD is released and activated by caspase cleavage of ICAD, and eventually causes degradation of nuclear DNA.</p> <p>Rabbit anti Murine CAD antibody recognizes a band of ~40 kDa in western blotting.</p>
Further Reading	<ol style="list-style-type: none"> 1. Enari, M. <i>et al.</i> (1998) A caspase-activated DNase that degrades DNA during apoptosis, and its inhibitor ICAD. Nature. 391 (6662): 43-50. 2. Sakahira, H. <i>et al.</i> (1998) Cleavage of CAD inhibitor in CAD activation and DNA degradation during apoptosis. Nature. 391 (6662): 96-9. 3. Liu, X. <i>et al.</i> (1997) DFF, a heterodimeric protein that functions downstream of caspase-3 to trigger DNA fragmentation during apoptosis. Cell. 89 (2): 175-84.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP486 10040</p>
Regulatory	For research purposes only

Related Products

[Recommended Secondary Antibodies](#)

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M364235:200529'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)