

Datasheet: AHP487

BATCH NUMBER 164759

Description:	RABBIT ANTI APAF1 (aa12-28)
Specificity:	APAF1 (aa12-28)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/1000 - 1/2000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species Human

Species Cross Reactivity Reacts with: Mouse
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to human Apaf-1 were raised by repeated immunisation of rabbits with highly

purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.02% Sodium Azide
---------------------------------	--------------------

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
---------------------------------------	-----------------------------

Immunogen	Synthetic peptide corresponding to amino acids 12-28 of human Apaf1.
------------------	--

External Database Links	UniProt: O14727 Related reagents Entrez Gene: 317 APAF1 Related reagents
--------------------------------	---

Synonyms	KIAA0413
-----------------	----------

RRID	AB_322893
-------------	-----------

Specificity	Rabbit anti Apaf1 antibody recognizes Apaf1 (Apoptosis protease activating factor 1), the mammalian homologue of the CED-4 cell death protein in <i>C. Elegans</i> . Apaf1 binds to cytochrome C and caspase-9 leading to caspase-9 activation, resulting in the beginning of the proteolytic cascade characteristic of apoptosis. Apaf1 can also associate with caspase-4 and caspase-8. Rabbit anti Apaf1 antibody (AHP487) recognizes a ~130 kDa band in western blots.
--------------------	--

Further Reading	1. Zou, H. <i>et al.</i> (1997) Apaf-1, a human protein homologous to <i>C. elegans</i> CED-4, participates in cytochrome c-dependent activation of caspase-3. Cell. 90 (3): 405-13. 2. Cecconi, F. <i>et al.</i> (1998) Apaf1 (CED-4 homolog) regulates programmed cell death in mammalian development. Cell. 94 (6): 727-37. 3. Hu, Y. <i>et al.</i> (1998) Bcl-XL interacts with Apaf-1 and inhibits Apaf-1-dependent caspase-9 activation. Proc Natl Acad Sci U S A. 95 (8): 4386-91.
------------------------	---

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP487 10040
--------------------------------------	---

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M382615:210513'

Printed on 12 Aug 2023
