

Datasheet: AHP627

BATCH NUMBER 165549

Description:	RABBIT ANTI HUMAN CD209
Specificity:	CD209
Other names:	DC-SIGN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1ug/ml - 2ug/ml
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Species Cross Reactivity	<p>Reacts with: Sheep</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
---------------------------------	---

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antisera to human DC-SIGN were raised by repeated immunisation of rabbits with highly

purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.02% Sodium Azide
---------------------------------	--------------------

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
---------------------------------------	-----------------------------

Immunogen	Synthetic peptide corresponding to amino acids 277-293 of human DC-SIGN.
------------------	--

External Database Links	UniProt: Q9NNX6 Related reagents Entrez Gene: 30835 CD209 Related reagents
--------------------------------	---

Synonyms	CLEC4L
-----------------	--------

RRID	AB_323683
-------------	-----------

Specificity	<p>Rabbit anti Human CD209 antibody recognizes an epitope within the extracellular domain of human DC-SIGN. DC-SIGN was designated CD209 at the 7th leucocyte typing workshop.</p> <p>Rabbit anti Human CD209 antibody detects a band of approximately 44 kDa by western blot.</p>
--------------------	---

References	<ol style="list-style-type: none">1. Akesson, C.P. <i>et al.</i> (2008) Phenotypic characterisation of intestinal dendritic cells in sheep. Dev Comp Immunol. 32: 837-49.2. Necchi, V. <i>et al.</i> (2009) Evidence for transepithelial dendritic cells in human <i>H. pylori</i> active gastritis. Helicobacter. 14: 208-22.3. Olivier, M. <i>et al.</i> (2012) Capacities of Migrating CD1b Lymph Dendritic Cells to Present <i>Salmonella</i> Antigens to Naive T Cells. PLoS One. 7: e30430.4. Emara, M. <i>et al.</i> (2012) Retagging identifies dendritic cell-specific intercellular adhesion molecule-3 (ICAM3)-grabbing non-integrin (DC-SIGN) protein as a novel receptor for a major allergen from house dust mite. J Biol Chem. 287 (8): 5756-63.5. Piercey Åkesson C <i>et al.</i> (2012) Phenotypic characterization of cells participating in transport of prion protein aggregates across the intestinal mucosa of sheep. Prion. 6 (3): 261-75.
-------------------	---

Further Reading	1. Geijtenbeek, T. B. <i>et al.</i> (2000) DC-SIGN, a dendritic cell-specific HIV-1-binding protein that enhances trans-infection of T cells. Cell. 100: 587-97.
------------------------	--

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP627 10040
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M382982:210513'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)