

Datasheet: AHP767

BATCH NUMBER 159671

Description:	RABBIT ANTI HUMAN EGF
Specificity:	EGF
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			5.0ug/ml
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			<0.1ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires heat-mediated pre-treatment of paraffin sections prior to staining. Citrate buffer pH 6.0 is recommended for this purpose

Target Species

Human

Product Form

Purified IgG - lyophilised

Reconstitution

Reconstitute with 0.1ml distilled water.

For long term storage the addition of 0.09% sodium azide is recommended.

N.B. For functional studies do not add azide.

Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.

Antiserum Preparation Antisera to human EGF were raised by repeated immunisations of rabbits with highly purified antigen.
Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers None present.

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0mg/ml after reconstitution

Immunogen Recombinant human EGF ([PHP030A](#)).

External Database Links

UniProt:

[P01133](#) [Related reagents](#)

Entrez Gene:

[1950](#) EGF [Related reagents](#)

RRID AB_2095966

Specificity **Rabbit anti Human EGF** recognizes human Epidermal Growth Factor, a potent stimulator and regulator of the proliferation of epidermal and epithelial cells, expressed by a variety of cells.

The effects of EGF are initiated through tyrosine kinase activity, following the binding of EGF to the extracellular domain of the ~170 kDa EGF receptor (EGFR), which also binds other EGF-like ligands with high affinity, including TGF- α (transforming growth factor alpha), VGF (vaccinia virus growth factor) and Betacellulin, expressed in the pancreas, small intestine and by certain tumour cells ([Seno *et al.* 1996](#)). Receptor activation through EGF binding triggers several signal transduction pathways, including the JAK/STAT and P13K/AKT pathways.

ELISA This purified human EGF antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with a [biotinylated human EGF antibody](#) (AHP767B) as the detection reagent and [recombinant human EGF](#) (PHP030A) as the standard.

Histology Positive Control Tissue Human kidney

Western Blotting Rabbit anti human EGF may be used in Western Blotting applications under either reducing or non-reducing conditions with [recombinant human EGF](#) (PHP030A) as the positive control.

Further Reading 1. Carpenter, G. & Cohen, S. (1990) Epidermal growth factor. [J Biol Chem. 265 \(14\)](#);

[7709-12.](#)

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10162 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP767>
10162

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M372537:200706'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)