

Datasheet: AHP831

BATCH NUMBER 161213

Description:	GOAT ANTI HUMAN VEGF
Specificity:	VEGF
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			5.0 - 15.0ug/ml
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.05 - 0.10ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)Heat-mediated antigen retrieval is required prior to staining paraffin sections with this antibody. Citrate buffer, pH6.0, is recommended for this purpose.

Target Species

Human

Product Form

Purified IgG - lyophilised

Reconstitution

Reconstitute with 0.1ml distilled water. For long term storage, the addition of 0.09% sodium azide is recommended. NB: For functional studies do not add azide.

Antiserum Preparation

Antisera to human VEGF were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution.
Immunogen	Recombinant human VEGF
External Database Links	<p>UniProt: P15692 Related reagents</p> <p>Entrez Gene: 7422 VEGFA Related reagents</p>
Synonyms	VEGF
RRID	AB_2212845
Specificity	<p>Goat anti Human VEGF antibody recognizes human VEGF (vascular endothelial growth factor), a heparin-binding cytokine and member of the PDGF/VEGF growth factor family, which exists as several splice variants of the VEGF gene. VEGF, originally named vascular permeability factor due to its ability to act as vasodilator and increase microvascular permeability, is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis.</p> <p>The functional forms of VEGF exist as disulphide-linked homodimers and are classified according to the number of amino acid residues within a constitutive polypeptide chain e.g VEGF₁₆₅. All VEGF family members act through binding to the tyrosine kinase receptors VEGFR-1 (Flt-1) and VEGFR-2 (KDR/CD309), whilst the closely-related VEGF-C and VEGF-D, can also bind to VEGFR-3 (Flt-4).</p>
ELISA	This purified human VEGF antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our biotinylated human VEGF antibody (AHP831B) as the detection reagent and recombinant human VEGF (PHP091) as the standard.
Histology Positive Control Tissue	Human breast invasive ductal carcinoma
Western Blotting	This antibody may be used in Western Blotting under either reducing or non-reducing conditions with recombinant human VEGF (PHP091) as the positive control.
References	1. Hjelmgren, O. <i>et al.</i> (2016) Increased Vascularization in the Vulnerable Upstream Regions of Both Early and Advanced Human Carotid Atherosclerosis. PLoS One. 11 (12): e0166918.

2. Gamliel, M. *et al.* (2018) Trained Memory of Human Uterine NK Cells Enhances Their Function in Subsequent Pregnancies. [Immunity. 48 \(5\): 951-962.e5.](#)

Further Reading 1. Neufeld, G. *et al.* (1999) Vascular endothelial growth factor (VEGF) and its receptors. [FASEB J. 13: 9-22.](#)

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10162 available at: <https://www.bio-rad-antibodies.com/SDS/AHP831>
10162

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

Recommended Positive Controls

[RECOMBINANT HUMAN VEGF \(PHP091\)](#)

Recommended Useful Reagents

[GOAT ANTI HUMAN VEGF:Biotin \(AHP831B\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M372566:200706'

Printed on 16 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)