

Datasheet: AHP851Z

BATCH NUMBER 157959

Description:	RABBIT ANTI BOVINE INTERLEUKIN-1 BETA:Preservative Free
Specificity:	IL-1 BETA
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1ug/ml - 10ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Species Cross Reactivity	Does not react with:Human, Mouse, Rat
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to bovine IL-1 beta were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography on protein G.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine IL-1 beta.
External Database Links	<p>UniProt: P09428 Related reagents</p> <p>Entrez Gene: 281251 IL1B Related reagents</p>
RRID	AB_2124486
Specificity	<p>Rabbit anti Bovine Interleukin-1 beta antibody recognizes bovine interleukin-1 beta (IL-1 beta).</p> <p>IL-1 beta is principally expressed by macrophages, and plays an important immunomodulatory role in the bovine immune response.</p> <p>Rabbit anti Bovine Interleukin-1 beta antibody recognizes both native and recombinant bovine IL-1 beta. No detectable cross-reactivity was observed with human, mouse or rat IL-1 beta.</p>
ELISA	This antibody may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA system, in combination with biotinylated Rabbit anti Bovine interleukin-21 antibody (AHP851B) as the detection antibody. Recombinant bovine interleukin-21 (PBP008) may be used as a standard.
References	<ol style="list-style-type: none"> 1. Wenz, J.R. <i>et al.</i> (2010) Factors associated with concentrations of select cytokine and acute phase proteins in dairy cows with naturally occurring clinical mastitis. J Dairy Sci. 93: 2458-70. 2. Rinaldi, M. <i>et al</i> (2010) A sentinel function for teat tissues in dairy cows: dominant innate immune response elements define early response to <i>E. coli</i> mastitis. Funct Integr Genomics. 10: 21-38. 3. Bougarn, S. <i>et al.</i> (2010) Muramyl dipeptide synergizes with <i>Staphylococcus aureus</i> lipoteichoic acid to recruit neutrophils in the mammary gland and to stimulate mammary epithelial cells. Clin Vaccine Immunol. 17 (11): 1797-809. 4. Simojoki, H. <i>et al.</i> (2011) Innate immune response in experimentally induced bovine intramammary infection with <i>Staphylococcus simulans</i> and <i>S. epidermidis</i>. Vet Res. 42: 49. 5. Harte,C. <i>et al.</i> (2017) Alum Activates the Bovine NLRP3 Inflammasome Frontiers in Immunology. 8 [Epub ahead of print].
Storage	<p>Store at -20°C only.</p> <p>This product should be stored undiluted.</p>

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/AHP851Z 10162
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[RABBIT ANTI BOVINE INTERLEUKIN-1 BETA:Biotin \(AHP851B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M364334:200529'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)