

Datasheet: AHP852Z

BATCH NUMBER 159804

Description:	RABBIT ANTI BOVINE TNF ALPHA:Preservative Free
Specificity:	TNF ALPHA
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1ug/ml - 10ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			1ug/ml - 10ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Bovine

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to bovine TNF alpha were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography on protein G.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers None present

Endotoxin Level <0.015 EU/ug

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine TNF alpha.
External Database Links	<p>UniProt: Q06599 Related reagents</p> <p>Entrez Gene: 280943 TNF Related reagents</p>
Synonyms	TNFA, TNFSF2
RRID	AB_567236
Specificity	<p>Rabbit anti Bovine TNF alpha antibody recognizes bovine tumour necrosis factor alpha (TNF alpha).</p> <p>TNF alpha is produced by a range of cells including macrophages, monocytes, T and B cells, neutrophils and endothelial cells. TNF alpha has a key regulatory role in the bovine immune response, where it is involved in modulating the expression of numerous cytokines, growth factors, inflammatory mediators and acute phase proteins.</p>
References	<ol style="list-style-type: none"> 1. Baravalle ME <i>et al.</i> (2015) Altered Expression of Pro-inflammatory Cytokines in Ovarian Follicles of Cows with Cystic Ovarian Disease. J Comp Pathol. 153 (2-3): 116-30. 2. Cigliano, L. <i>et al.</i> (2016) Evaluation of serum markers of blood redox homeostasis and inflammation in PCB naturally contaminated heifers undergoing decontamination. Sci Total Environ. 542 (Pt A): 653-64. 3. Canal AM <i>et al.</i> (2017) Immunohistochemical detection of pro-inflammatory and anti-inflammatory cytokines in granulomas in cattle with natural <i>Mycobacterium bovis</i> infection. Res Vet Sci. 110: 34-39. 4. Castel, A. <i>et al.</i> (2021) Recording and manipulation of vagus nerve electrical activity in chronically instrumented unanesthetized near term fetal sheep. J Neurosci Methods. 360: 109257.
Storage	<p>Store at -20°C only.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/AHP852Z</p> <p>10162</p>

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M364336:200529'

Printed on 19 Jan 2024