

Datasheet: AHP865

BATCH NUMBER 157244

Description:	GOAT ANTI PIG IgG (H/L)
Specificity:	IgG (H/L)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/100 - 1/10,000
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Pig

Product Form

Purified IgG - liquid

Antiserum Preparation

Antisera to porcine IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography using antigen coupled to agarose beads.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Approx. Protein

IgG concentration 1.0 mg/ml

Concentrations

RRID AB_2121737

Specificity **Goat anti Pig IgG (H/L) antibody** recognizes porcine IgG (heavy and light chains) and has been shown to react specifically with porcine IgG and with light chains common to other pig immunoglobulins, by immunoelectrophoresis and ELISA.

Goat anti Pig IgG (H/L) antibody may cross react with IgG from other species.

References

1. Tian, F. *et al.* (2010) Immune Events Associated with High Level Protection against *Schistosoma japonicum* Infection in Pigs Immunized with UV-Attenuated Cercariae. [PLoS One. 5\(10\): e13408.](#)
2. Lin, D. *et al.* (2011) Multiple vaccinations with UV- attenuated cercariae in pig enhance protective immunity against *Schistosoma japonicum* infection as compared to single vaccination. [Parasit Vectors. 4:103.](#)
3. Olasz, F. *et al.* (2016) Immunological and biochemical characterisation of 7ap, a short protein translated from an alternative frame of ORF7 of PRRSV. [Acta Vet Hung. 64 \(2\): 273-87.](#)
4. Rungelrath, V. *et al.* (2018) IgM cleavage by *Streptococcus suis*. reduces IgM bound to the bacterial surface and is a novel complement evasion mechanism. [Virulence. 9 \(1\): 1314-1337.](#)
5. Grzeškowiak, Ł. *et al.* (2022) *Clostridioides difficile*-mesocolonic oedema in neonatal suckling piglets develops regardless of the fibre composition in sow's diets. [Animal. 17 \(2\): 100697.](#)

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP865>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M364346:200529'

