

Datasheet: AHP865P

| | |
|----------------------|-----------------------------|
| Description: | GOAT ANTI PIG IgG (H/L):HRP |
| Specificity: | IgG (H/L) |
| Format: | HRP |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|----------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/10,000 - 1/100,000 |
| Immunoprecipitation | | | ▪ | |
| Western Blotting | | | ▪ | |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|---|
| Target Species | Pig |
| Product Form | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid |
| Antiserum Preparation | Antisera to porcine IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography. |
| Buffer Solution | Phosphate buffered saline. |
| Preservative | 0.05% ProClin 300 |
| Stabilisers | 0.2% Bovine Serum Albumin |
| Approx. Protein Concentrations | IgG concentration 1.0 mg/ml. |

RRID AB_2121735.

Specificity **Goat anti Pig IgG (H/L) antibody** recognizes porcine IgG (heavy and light chains) and has been shown to react specifically with porcine IgG and with light chains common to other pig immunoglobulins, by immunoelectrophoresis and ELISA.

Goat anti Pig IgG (H/L) antibody may cross react with IgG from other species.

- References**
1. Tian, F. *et al.* (2010) Immune Events Associated with High Level Protection against *Schistosoma japonicum* Infection in Pigs Immunized with UV-Attenuated Cercariae. [PLoS One. 5\(10\): e13408.](#)
 2. Lin, D. *et al.* (2011) Multiple vaccinations with UV- attenuated cercariae in pig enhance protective immunity against *Schistosoma japonicum* infection as compared to single vaccination. [Parasit Vectors. 4:103.](#)
 3. Olasz, F. *et al.* (2016) Immunological and biochemical characterisation of 7ap, a short protein translated from an alternative frame of ORF7 of PRRSV. [Acta Vet Hung. 64 \(2\): 273-87.](#)
 4. Rungelrath, V. *et al.* (2018) IgM cleavage by *Streptococcus suis*. reduces IgM bound to the bacterial surface and is a novel complement evasion mechanism. [Virulence. 9 \(1\): 1314-1337.](#)
 5. Grześkowiak, Ł. *et al.* (2022) *Clostridioides difficile*-mesocolonic oedema in neonatal suckling piglets develops regardless of the fibre composition in sow's diets. [Animal. 17 \(2\): 100697.](#)
-

Storage Store at +4C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch.

Acknowledgements ProClin is a trademark of The Dow Chemical Company ("Dow") or an affiliated company of Dow.

Health And Safety Information Material Safety Datasheet documentation #20391 available at: <https://www.bio-rad-antibodies.com/SDS/AHP865P>
20391

Regulatory For research purposes only.

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M432723:240912'

Printed on 12 Sep 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)