

Datasheet: AHP932

BATCH NUMBER 164895

Description:	RABBIT ANTI HUMAN DOPAMINE TRANSPORTER (C-TERMINAL)
Specificity:	DOPAMINE TRANSPORTER (C-TERMINAL)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/1000
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/1000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Monkey, Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antisera to human dopamine transporter (CT) were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	10mM HEPES pH7.5 150mM NaCl
Preservative Stabilisers	0.09% Sodium Azide 0.01% Bovine Serum Albumin 50% Glycerol
Immunogen	Keyhole limpet haemocyanin conjugated synthetic peptide corresponding to the C-terminus region (CT) of human dopamine transporter.
External Database Links	UniProt: Q01959 Related reagents Entrez Gene: 6531 SLC6A3 Related reagents
Synonyms	DAT1
RRID	AB_2286179
Specificity	Rabbit anti Human dopamine transporter antibody recognizes an epitope within the C-terminal region of dopamine transporter (DAT), also known as solute carrier family 6. DAT is principally expressed in the presynaptic membrane of dopaminergic nerve terminals, where it mediates the rapid removal of dopamine from the synaptic cleft after neuronal stimulation. As such, DAT plays an important role in regulating the availability of dopamine in the synaptic cleft and thus in modulating the physiological functions of dopamine. The reuptake of dopamine by DAT is inhibited by psychostimulants such as cocaine and amphetamine.
Western Blotting	AHP932 detects a band of approximately 88kDa in human caudate cell lysates.
References	1. Mozley, L.H. <i>et al.</i> (2001) Striatal dopamine transporters and cognitive functioning in healthy men and women. Am J Psychiatry. 158 (9): 1492-9.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety	Material Safety Datasheet documentation #10088 available at:

Information <https://www.bio-rad-antibodies.com/SDS/AHP932>
10088

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M382995:210513'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)