

Datasheet: AHP945Z

BATCH NUMBER 157497

Description:	RABBIT ANTI CHICKEN INTERFERON GAMMA:Preservative Free
Specificity:	IFN GAMMA
Other names:	INTERFERON GAMMA
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			
Functional Assays	▪			10ug/ml - 50ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Chicken
-----------------------	---------

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antisera to chicken Interferon Gamma were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography on protein A.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	None present.
---------------------------------	---------------

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant chicken interferon gamma
External Database Links	<p>UniProt: P49708 Related reagents</p> <p>Entrez Gene: 396054 IFNG Related reagents</p>
RRID	AB_566996
Specificity	<p>Rabbit anti Chicken interferon gamma antibody recognizes chicken interferon gamma (chIFN-gamma). There are at least three distinct chicken interferon proteins, designated as chIFN-alpha, chIFN-beta and chIFN-gamma.</p> <p>chIFN-gamma was first identified, when supernatant from chicken spleen cells stimulated with concanavalin A (ConA), displayed macrophage-activating factor (MAF) activity - the principle function of mammalian type II IFN-gamma (Vonbulow <i>et al.</i> 1984). Cloning of IFN-gamma cDNA from activated chicken T cells indicates that recombinant chicken IFN-gamma activates macrophages, strongly upregulating the expression of MHC class II and iNOS genes (Weining <i>et al.</i> 1996).</p> <p>Rabbit anti Chicken interferon gamma antibody can be used to neutralize chIFN-gamma activity and to discriminate the antiviral activity of chIFN-gamma from that of chIFN-alpha and chIFN-beta (Weining <i>et al.</i> 1996).</p>
Western Blotting	AHP945Z detects a band of approximately 18 kDa.
References	<ol style="list-style-type: none"> Wei, Y. <i>et al.</i> (2016) Generation and protective efficacy of a cold-adapted attenuated avian H9N2 influenza vaccine. Sci Rep. 6: 30382. Andersen, S.H. <i>et al.</i> (2017) Quantification and phenotypic characterisation of peripheral IFN-γ producing leucocytes in chickens vaccinated against Newcastle disease. Vet Immunol Immunopathol. 193-194: 18-28. Wang, M. <i>et al.</i> (2018) Cross- immunity of a H9N2 live attenuated influenza vaccine against H5N2 highly pathogenic avian influenza virus in chickens. Vet Microbiol. 220: 57-66.
Further Reading	<ol style="list-style-type: none"> Schultz, U. <i>et al.</i> (1995) Recombinant chicken interferon: a potent antiviral agent that lacks intrinsic macrophage activating factor activity. Eur J Immunol. 25 (3): 847-51. Weining, K.C. <i>et al.</i> (1996) Biological properties of recombinant chicken interferon-gamma. Eur J Immunol. 26 (10): 2440-7. Staeheli, P. <i>et al.</i> (2001) Cytokines of birds: conserved functions--a largely different look. J Interferon Cytokine Res. 21 (12): 993-1010.
Storage	Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/AHP945Z 10162
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M386898:210601'

Printed on 19 Jan 2024