

Datasheet: AHP948P

BATCH NUMBER 158415

| | |
|----------------------|-------------------------|
| Description: | SHEEP ANTI DOG IgG2:HRP |
| Specificity: | IgG2 |
| Format: | HRP |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|----------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | ▪ | | | 1/200 - 1/500 |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/10,000 - 1/100,000 |
| Immunoprecipitation | | | ▪ | |
| Western Blotting | ▪ | | | 1/1000 - 1/10,000 |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Dog

Product Form

Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation

Antisera to canine IgG2 were raised by repeated immunisations of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.05% Proclin™ 300
0.2% Bovine Serum Albumin

Approx. Protein Concentrations

IgG concentration 1.0 mg/ml

RRID AB_567050

Specificity **Sheep anti Dog IgG2 polyclonal antibody** recognizes canine IgG2 by electrophoresis and ELISA. No cross reactivity was detected against other canine immunoglobulin classes, isolated immunoglobulin light chains or no immunoglobulin proteins.

Goat anti Dog IgG2 may cross react with IgG2 from other species.

References

1. Bird, R.C. *et al.* (2011) An autologous dendritic cell canine mammary tumor hybrid-cell fusion vaccine. [Cancer Immunol Immunother. 60 \(1\): 87-97.](#)
2. Agallou M *et al.* (2016) Identification of Immunoreactive *Leishmania infantum* Protein Antigens to Asymptomatic Dog Sera through Combined Immunoproteomics and Bioinformatics Analysis. [PLoS One. 11 \(2\): e0149894.](#)
3. Martínez Abad, L.P. *et al.* (2017) Diagnostic accuracy of rKLO8 versus rK26 ELISAs for screening of canine visceral leishmaniasis. [Acta Trop. 166: 133-8.](#)
4. Khantavee, N. *et al.* (2020) Antibody levels to *Malassezia pachydermatis*. and *Staphylococcus pseudintermedius*. in atopic dogs and their relationship with lesion scores. [Vet Dermatol. 31 \(2\): 111-115.](#)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Acknowledgements Proclin™ 300 is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow.

Health And Safety Information Material Safety Datasheet documentation #20391 available at: <https://www.bio-rad-antibodies.com/SDS/AHP948P>
20391

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M364385:200529'

