

Datasheet: AHP949P

BATCH NUMBER 152027

Description:	RABBIT ANTI SHEEP IgA:HRP
Specificity:	IgA
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/10,000 - 1/100,000
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Sheep
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Antiserum Preparation	Antisera to ovine IgA were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	0.05% Proclin™ 300
Stabilisers	0.2% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

RRID	AB_567005
Specificity	<p>Rabbit anti sheep IgA antibody recognizes ovine IgA. This antibody has been shown to react with ovine IgA by immunoelectrophoresis and ELISA.</p> <p>Rabbit anti sheep IgA antibody may cross react with IgA from other species.</p>
References	<ol style="list-style-type: none"> Murphy, L. <i>et al.</i> (2010) Genetic variation among lambs in peripheral IgE activity against the larval stages of <i>Teladorsagia circumcincta</i>. Parasitology. 137: 1249-60. Hine, B.C. <i>et al.</i> (2010) Selective transport of IgE into ovine mammary secretions. Res Vet Sci. 89 (2): 184-90. Vande Walle, K. <i>et al.</i> (2011) Rectal inoculation of sheep with <i>E. coli</i> O157:H7 results in persistent infection in the absence of a protective immune response. Vet Microbiol. 147 (3-4): 376-82. Díaz, A.G. <i>et al.</i> (2016) Spray dried microspheres based on chitosan: A promising new carrier for intranasal administration of polymeric antigen BLSOmp31 for prevention of ovine brucellosis. Mater Sci Eng C Mater Biol Appl. 62: 489-96. Díaz, A.G. <i>et al.</i> (2016) Immune response induced by conjunctival immunization with polymeric antigen BLSOmp31 using a thermoresponsive and mucoadhesive in situ gel as vaccine delivery system for prevention of ovine brucellosis. Vet Immunol Immunopathol. 178: 50-6. Subharat, S. <i>et al.</i> (2016) Vaccination of Sheep with a Methanogen Protein Provides Insight into Levels of Antibody in Saliva Needed to Target Ruminal Methanogens. PLoS One. 11 (7): e0159861. Watt, K.A. <i>et al.</i> (2016) Fecal antibody levels as a noninvasive method for measuring immunity to gastrointestinal nematodes in ecological studies. Ecol Evol. 6 (1): 56-67. Sparks, A.M. <i>et al.</i> (2018) Natural Selection on Antihelminth Antibodies in a Wild Mammal Population. Am Nat. 192 (6): 745-60. Sparks, A.M. <i>et al.</i> (2019) The genetic architecture of helminth-specific immune responses in a wild population of Soay sheep (<i>Ovis aries</i>). PLoS Genet. 15 (11): e1008461. Ramos, A. <i>et al.</i> (2018) Melatonin enhances responsiveness to Dichelobacter nodosus vaccine in sheep and increases peripheral blood CD4 T lymphocytes and IgG-expressing B lymphocytes. Vet Immunol Immunopathol. 206: 1-8. Fairlie-Clarke, K. <i>et al.</i> (2019) Salivary IgA: A biomarker for resistance to <i>Teladorsagia circumcincta</i> and a new estimated breeding value. Vet Parasitol. 269: 16-20.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Acknowledgements	Proclin™ 300 is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow.
Health And Safety Information	<p>Material Safety Datasheet documentation #20391 available at: https://www.bio-rad-antibodies.com/SDS/AHP949P</p> <p>20391</p>

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M364387:200529'

Printed on 18 Jan 2024