

Datasheet: AHP950

BATCH NUMBER 171009

Description:	RABBIT ANTI SHEEP IgM
Specificity:	IgM
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/100 - 1/10,000
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Sheep

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to ovine IgM were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Specificity **Rabbit anti sheep IgM antibody** recognizes ovine immunoglobulin M.

The antibody has been shown to react specifically with ovine IgM by immunoelectrophoresis and ELISA.

The antibody may cross react with IgM from other species.

References

- Hine, B.C. *et al.* (2010) Selective transport of IgE into ovine mammary secretions. [Res Vet Sci. 89 \(2\): 184-90.](#)
- Vande Walle, K. *et al.* (2011) Rectal inoculation of sheep with *E. coli* O157:H7 results in persistent infection in the absence of a protective immune response. [Vet Microbiol. 147 \(3-4\): 376-82.](#)
- Marques, P.X. *et al.* (2011) Amniotic and allantoic fluids from experimentally infected sheep contain immunoglobulin specific for *Chlamydophila abortus*. [Vet Immunol Immunopathol. 140: 1-9.](#)
- Verhelst, D. *et al.* (2015) Interferon-gamma expression and infectivity of *Toxoplasma* infected tissues in experimentally infected sheep in comparison with pigs. [Vet Parasitol. 207 \(1-2\): 7-16.](#)
- Helder, M.R.K. *et al.* (2017) Xenoantigenicity of porcine decellularized valves. [J Cardiothorac Surg. 12 \(1\): 56.](#)

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP950>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M364388:200529'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)