

Datasheet: AHP950F

BATCH NUMBER 163650

Description:	RABBIT ANTI SHEEP IgM:FITC
Specificity:	IgM
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/200 - 1/2000
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Sheep		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

Antiserum Preparation Antisera to ovine IgM were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative 0.09% Sodium Azide
Stabilisers 0.2% Bovine Serum Albumin

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
RRID	AB_567060
Specificity	<p>Rabbit anti sheep IgM antibody recognizes ovine immunoglobulin M.</p> <p>The antibody has been shown to react specifically with ovine IgM by immunoelectrophoresis and ELISA.</p> <p>The antibody may cross react with IgM from other species.</p>
References	<ol style="list-style-type: none"> Hine, B.C. <i>et al.</i> (2010) Selective transport of IgE into ovine mammary secretions. Res Vet Sci. 89 (2): 184-90. Vande Walle, K. <i>et al.</i> (2011) Rectal inoculation of sheep with <i>E. coli</i> O157:H7 results in persistent infection in the absence of a protective immune response. Vet Microbiol. 147 (3-4): 376-82. Marques, P.X. <i>et al.</i> (2011) Amniotic and allantoic fluids from experimentally infected sheep contain immunoglobulin specific for <i>Chlamydophila abortus</i>. Vet Immunol Immunopathol. 140: 1-9. Verhelst, D. <i>et al.</i> (2015) Interferon-gamma expression and infectivity of <i>Toxoplasma</i> infected tissues in experimentally infected sheep in comparison with pigs. Vet Parasitol. 207 (1-2): 7-16. Helder, M.R.K. <i>et al.</i> (2017) Xenoantigenicity of porcine decellularized valves. J Cardiothorac Surg. 12 (1): 56.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/AHP950F</p> <p>10041</p>
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M364389:200529'

Printed on 18 Jan 2024