

Datasheet: AHP964

BATCH NUMBER 159754

Description:	RABBIT ANTI CASPASE-4 (N-TERMINAL)
Specificity:	CASPASE-4 (N-TERMINAL)
Other names:	ICH-2
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin	▪			2.0ug/ml
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Antisera to human Caspase-4 (NT) were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	A 16 amino acid peptide from human Caspase-4 amino terminus.
External Database Links	<p>UniProt: P49662 Related reagents</p> <p>Entrez Gene: 837 CASP4 Related reagents</p>
Synonyms	ICH2
RRID	AB_2070198
Specificity	<p>Rabbit anti caspase-4 (N-terminal) antibody recognizes an epitope within the N-terminal region (NT) of Caspase-4, otherwise known as ICH-2. Caspase-4, a member of the Caspase-1 subfamily of cysteine proteases, exists as an inactive pro-enzyme which undergoes proteolytic cleavage of its p30 precursor form, into smaller p20 and p10 subunits.</p> <p>Caspase-4 is localized to the endoplasmic reticulum (ER) membrane and is involved in response to ER membrane stress in conditions including Alzheimer's disease and Ischemia and confirms the involvement of Caspase-4 as an instigator of cellular apoptosis.</p>
Immunohistology	*This product does not require protein digestion pre-treatment of paraffin embedded sections e.g. trypsin or pronase prior to staining. This product does not require antigen retrieval using heat treatment prior to staining of paraffin sections.
Histology Positive Control Tissue	Spleen tissue.
Western Blotting	AHP964 detects a band of approximately 34 kDa in Ramos cell lysates.
References	1. Katayama, T. <i>et al.</i> (2004) Induction of neuronal death by ER stress in Alzheimer's disease. J Chem Neuroanat. 28 (1-2): 67-78.
Further Reading	<p>1. Kamens, J. <i>et al.</i> (1995) Identification and characterization of ICH-2, a novel member of the interleukin-1 beta-converting enzyme family of cysteine proteases. J Biol Chem. 270 (25): 15250-6.</p> <p>2. Cohen, G.M. (1997) Caspases: the executioners of apoptosis. Biochem J. 326 (Pt 1): 1-16.</p>

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP964 10040
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

- [ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
- [TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M382999:210513'

Printed on 18 Jan 2024