

Datasheet: AHP966

BATCH NUMBER 158967

Description:	RABBIT ANTI CASPASE-7 (C-TERMINAL)
Specificity:	CASPASE-7 (C-TERMINAL)
Other names:	ICE-LAP3
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			2.0ug/ml
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1.0ug/ml - 2.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Mouse, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid

Antiserum Preparation Antisera to human Caspase-7 (CT) were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.02% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Immunogen A 16 amino acid peptide from human Caspase-7 carboxy-terminus

External Database Links

UniProt:
[P55210](#) [Related reagents](#)

Entrez Gene:
[840](#) CASP7 [Related reagents](#)

Synonyms MCH3

RRID AB_566403

Specificity **Rabbit anti Human caspase-7 antibody** recognizes an epitope within the C-terminal region (CT) of Caspase-7, a ~35 kDa cysteine protease, otherwise known as ICE-like Apoptotic Protease 3 (ICE-LAP3).

Caspase-7, a member of the ICE/Ced-3 subfamily, is an executioner caspase which undergoes proteolytic cleavage of its precursor to form active p12 and p20 subunits. Evidence that activation of Caspase-7 occurs during cell death induced by the cytokine death receptors Fas/APO-1 and the receptor of tumour necrosis factor (TNFR-1), coupled with the fact that granzyme-B activated Caspase-7 cleaves the nuclear enzyme poly (ADP-ribose) polymerase (PARP), suggests an important role for Caspase-7 in both cytokine-mediated and granzyme-B mediated apoptosis pathways.

Histology Positive Control Tissue Skeletal muscle tissue

Further Reading

1. Duan, H. *et al.* (1996) ICE-LAP3, a novel mammalian homologue of the *Caenorhabditis elegans* cell death protein Ced-3 is activated during Fas- and tumor necrosis factor-induced apoptosis. [J Biol Chem. 271 \(3\): 1621-5.](#)
2. Juan, T.S. *et al.* (1997) Identification and mapping of Casp7, a cysteine protease resembling CPP32 beta, interleukin-1 beta converting enzyme, and CED-3. [Genomics. 40 \(1\): 86-93.](#)
3. Cohen, G.M. (1997) Caspases: the executioners of apoptosis. [Biochem J. 326 \(Pt 1\): 1-16.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at

-20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP966 10040
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

- [ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
- [TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M383000:210513'

Printed on 18 Jan 2024