

Datasheet: AHP966

**BATCH NUMBER 170214**

<b>Description:</b>	RABBIT ANTI CASPASE-7 (C-TERMINAL)
<b>Specificity:</b>	CASPASE-7 (C-TERMINAL)
<b>Other names:</b>	ICE-LAP3
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			2.0ug/ml
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1.0ug/ml - 2.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	<p>Reacts with: Mouse, Rat</p> <p><b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
<b>Product Form</b>	Purified IgG - liquid

**Antiserum Preparation** Antisera to human Caspase-7 (CT) were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

---

**Buffer Solution** Phosphate buffered saline

---

**Preservative Stabilisers** 0.02% Sodium Azide

---

**Approx. Protein Concentrations** IgG concentration 1.0 mg/ml

---

**Immunogen** A 16 amino acid peptide from human Caspase-7 carboxy-terminus

---

**External Database Links**

**UniProt:**  
[P55210](#)    [Related reagents](#)

**Entrez Gene:**  
[840](#)    CASP7    [Related reagents](#)

---

**Synonyms** MCH3

---

**RRID** AB\_566403

---

**Specificity** **Rabbit anti Human caspase-7 antibody** recognizes an epitope within the C-terminal region (CT) of Caspase-7, a ~35 kDa cysteine protease, otherwise known as ICE-like Apoptotic Protease 3 (ICE-LAP3).

Caspase-7, a member of the ICE/Ced-3 subfamily, is an executioner caspase which undergoes proteolytic cleavage of its precursor to form active p12 and p20 subunits. Evidence that activation of Caspase-7 occurs during cell death induced by the cytokine death receptors Fas/APO-1 and the receptor of tumour necrosis factor (TNFR-1), coupled with the fact that granzyme-B activated Caspase-7 cleaves the nuclear enzyme poly (ADP-ribose) polymerase (PARP), suggests an important role for Caspase-7 in both cytokine-mediated and granzyme-B mediated apoptosis pathways.

---

**Histology Positive Control Tissue** Skeletal muscle tissue

---

**Further Reading**

1. Duan, H. *et al.* (1996) ICE-LAP3, a novel mammalian homologue of the *Caenorhabditis elegans* cell death protein Ced-3 is activated during Fas- and tumor necrosis factor-induced apoptosis. [J Biol Chem. 271 \(3\): 1621-5.](#)
2. Juan, T.S. *et al.* (1997) Identification and mapping of Casp7, a cysteine protease resembling CPP32 beta, interleukin-1 beta converting enzyme, and CED-3. [Genomics. 40 \(1\): 86-93.](#)
3. Cohen, G.M. (1997) Caspases: the executioners of apoptosis. [Biochem J. 326 \( Pt 1\): 1-16.](#)

---

**Storage** Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP966">https://www.bio-rad-antibodies.com/SDS/AHP966</a> 10040
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

### Recommended Useful Reagents

- [ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
- [TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M364402:200529'

**Printed on 25 Mar 2023**