

Datasheet: BLP086

**BATCH NUMBER 163754**

<b>Description:</b>	SpyTag3 PEPTIDE
<b>Name:</b>	SPYTAG3 PEPTIDE
<b>Format:</b>	Purified
<b>Product Type:</b>	Blocking Peptide
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Protein/peptide tag
<b>Product Form</b>	Purified peptide - liquid
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	1 mM (1.93 mg/ml)

**Product Information** **The SpyTag3 peptide** is provided as a blocking peptide for use with antibodies ligated to SpyCatcher, BiCatcher and FcCatcher products. SpyTag3 peptide binds to excess SpyCatcher sites that have not reacted with an antibody during ligation. The blocking peptide will bind to SpyCatcher versions 1, 2 and 3.

A recombinant antibody or any other protein with a SpyTag can be ligated to one of our many prefabricated SpyCatcher protein modules. The reaction between SpyTag and SpyCatcher requires only mixing and a short incubation. If the ligation reaction mixture

contains an excess of the SpyCatcher partner, unreacted SpyCatcher sites remaining in the reaction mixture could interfere with other assay components and generate background signal. This can be prevented by the addition of SpyTag3 peptide to the final coupled antibody or protein preparation to block unreacted SpyCatcher sites.

It is recommended to use a 2.5-fold molar excess of SpyTag3 peptide over SpyCatcher sites, i.e. 2.5-fold excess for SpyCatcher-coupled antibodies and 5-fold excess for BiCatcher- and FcCatcher-coupled antibodies and incubate for at least 5 minutes.

Visit [bio-rad-antibodies.com/SpyCatcher](http://bio-rad-antibodies.com/SpyCatcher) for further information.

---

<b>Purity</b>	>90% by RP-HPLC
---------------	-----------------

---

<b>Storage</b>	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Acknowledgements</b>	Sold under license of patent and patent applications, including but not limited to the international PCT application PCT/GB2018/051065, the European patent EP2534484B1, UK patent GB1706430.4, the US patents US9,547,003, US10,247,727, US patent application US16/274908 and corresponding patents
-------------------------	---

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10209 available at: <a href="https://www.bio-rad-antibodies.com/SDS/BLP086">https://www.bio-rad-antibodies.com/SDS/BLP086</a> 10209
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Useful Reagents

[SpyCatcher2 \(TZC001\)](#)

[SpyCatcher2:Biotin \(TZC001B\)](#)

[SpyCatcher2-CYS \(TZC001CYS\)](#)

[BiSpyCatcher2 \(TZC002\)](#)

[BiSpyCatcher2:Biotin \(TZC002B\)](#)

[BiSpyCatcher2-CYS \(TZC002CYS\)](#)

[BiSpyCatcher2-CYS3 \(TZC002CYS3\)](#)

[BiSpyCatcher2:HRP \(TZC002P\)](#)

[SpyCatcher3 \(TZC025\)](#)

[SpyCatcher3-CYS \(TZC025CYS\)](#)

[HUMAN IgG1-FcSpyCatcher3 \(TZC009\)](#)

[MOUSE IgG2a-FcSpyCatcher3 \(TZC012\)](#)

[RABBIT IgG-FcSpyCatcher3 \(TZC013\)](#)

[HUMAN IgG2-FcSpyCatcher3 \(TZC016\)](#)

[HUMAN IgG3-FcSpyCatcher3 \(TZC017\)](#)

[HUMAN IgG4-FcSpyCatcher3 \(TZC018\)](#)

[HUMAN IgG4-Pro-FcSpyCatcher3 \(TZC019\)](#)

[HUMAN IgA-FcSpyCatcher3 \(TZC020\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M399023:220627'

**Printed on 26 Jun 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)