



Datasheet: BUF017B

BATCH NUMBER 171813

Description:	PEROXIDE BLOCKING REAGENT
Name:	PEROXIDE BLOCKING REAGENT
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	50 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Intended Use

BUF017B is a 3% (v/v) solution of hydrogen peroxide. This solution may be used to block endogenous peroxidase activity in paraffin embedded tissue sections. Incubate re-hydrated sections in this buffer for 15 minutes, and rinse well with water before continuing with staining protocol.

Storage

This product is shipped at ambient temperature.
Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10233 available at: <https://www.bio-rad-antibodies.com/SDS/BUF017B>

Regulatory

For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M442605:250528'

