

## Datasheet: BUF030C

**BATCH NUMBER 162555**

<b>Description:</b>	5x ELISA COATING BUFFER
<b>Name:</b>	5x ELISA COATING BUFFER
<b>Format:</b>	5 X Concentrate
<b>Product Type:</b>	Accessory Reagent
<b>Quantity:</b>	1000 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only.

<b>Product Form</b>	5 x concentrate - liquid
---------------------	--------------------------

<b>Preservative Stabilisers</b>	0.5% Sodium Azide
---------------------------------	-------------------

<b>Intended Use</b>	<p>BUF030C is a proprietary pH 7.8 - 8.2 coating buffer for use in all ELISA formats. Designed to stabilize the tertiary structure of the antibody or antigen and maximize adsorption onto polystyrene plates. It is especially effective in sandwich ELISA, preserving the antigen recognition regions of the coating antibody, allowing greater binding reactivity.</p>
---------------------	---

BUF030C contains an antimicrobial agent allowing coating to be performed at room temperature. Dried plates can be kept at 4°C for long term storage.

- |                             |  |
|-----------------------------|--|
| <b>Instructions For Use</b> | <ol style="list-style-type: none"><li>1. Dilute the required volume of BUF030C 1:5 with deionised water and mix for 5-15 minutes to prepare the working solution.</li><li>2. Add antibody or antigen to the working solution and mix for 5-15 minutes.</li><li>3. Coat plate (suggested volume 100-200µl per well) and incubate at room temperature for 6-24 hours, protected from light and covered to avoid evaporation.</li></ol> |
|-----------------------------|--|

**References**

1. Böröcz, K. *et al.* (2022) Dynamic Features of Herd Immunity: Similarities in Age-Specific Anti-Measles Seroprevalence Data between Two Countries of Different Epidemiological History. [J Clin Med. 11 \(4\): 1145.](#)
2. Hayden, Z. *et al.* (2021) Toll-Like Receptor Homolog CD180 Expression Is Diminished on Natural Autoantibody-Producing B Cells of Patients with Autoimmune CNS Disorders. [J Immunol Res. 2021: 9953317.](#)

**Storage**

Store at +4°C. DO NOT FREEZE.  
 This product should be stored undiluted. Should this product contain a precipitate, heat gently until dissolved. Do not boil.

**Guarantee**

Guaranteed until date of expiry. Please see product label.

**Health And Safety Information**

Material Safety Datasheet documentation #10383 available at:  
 10383: <https://www.bio-rad-antibodies.com/uploads/MSDS/10383.pdf>

**Regulatory**

For research purposes only

## Related Products

### Recommended Useful Reagents

[10x ELISA WASH BUFFER \(BUF031A\)](#)

[ELISA BSA BLOCK \(BUF032A\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
 'M395203:220323'

**Printed on 08 Feb 2023**