



## Datasheet: BUF036A

**BATCH NUMBER 161714**

|                      |                               |
|----------------------|-------------------------------|
| <b>Description:</b>  | 10X PHOSPHATE BUFFERED SALINE |
| <b>Name:</b>         | 10X PHOSPHATE BUFFERED SALINE |
| <b>Format:</b>       | 10 X Concentrate              |
| <b>Product Type:</b> | Accessory Reagent             |
| <b>Quantity:</b>     | 1 LITRE                       |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|             | Yes | No | Not Determined | Suggested Dilution |
|-------------|-----|----|----------------|--------------------|
| ELISA       | ▪   |    |                |                    |
| Immunoassay | ▪   |    |                | 1/10               |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

|                                      |   |
|--------------------------------------|---|
| <b>Product Form</b>                  | 10 x concentrate - liquid pH7.3 +/-0.1  |
| <b>Intended Use</b>                  | Dilute 1:9 with distilled water to prepare solution of 0.01M PBS for use in washing of sections in immunohistology staining procedures.   |
| <b>Storage</b>                       | Store at room temperature   |
| <b>Guarantee</b>                     | Guaranteed until date of expiry. Please see product label.  |
| <b>Health And Safety Information</b> | Material Safety Datasheet documentation #10216 available at: <a href="https://www.bio-rad-antibodies.com/SDS/BUF036A">https://www.bio-rad-antibodies.com/SDS/BUF036A</a><br>10216 |
| <b>Regulatory</b>                    | For research purposes only  |

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M390133:210825'

**Printed on 12 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)