



Datasheet: BUF036A

BATCH NUMBER 172483

| | |
|----------------------|-------------------------------|
| Description: | 10X PHOSPHATE BUFFERED SALINE |
| Name: | 10X PHOSPHATE BUFFERED SALINE |
| Format: | 10 X Concentrate |
| Product Type: | Accessory Reagent |
| Quantity: | 1 LITRE |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Immunoassay | ▪ | | | 1/10 |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|--------------------------------------|--|
| Product Form | 10 x concentrate - liquid pH7.3 +/-0.1 |
| Intended Use | Dilute 1:9 with distilled water to prepare solution of 0.01M PBS for use in washing of sections in immunohistology staining procedures. |
| Storage | This product is shipped at ambient temperature. Store at room temperature |
| Guarantee | Guaranteed until date of expiry. Please see product label. |
| Health And Safety Information | Material Safety Datasheet documentation #10216 available at: https://www.bio-rad-antibodies.com/SDS/BUF036A |
| Regulatory | For research purposes only |

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M442609:250528'

