



Datasheet: BUF036A

BATCH NUMBER 172511

Description:	10X PHOSPHATE BUFFERED SALINE
Name:	10X PHOSPHATE BUFFERED SALINE
Format:	10 X Concentrate
Product Type:	Accessory Reagent
Quantity:	1 LITRE

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	10 x concentrate - liquid pH7.3 +/-0.1
Intended Use	Dilute 1:9 with distilled water to prepare solution of 0.01M PBS for use in washing of sections in immunohistology staining procedures.
Storage	This product is shipped at ambient temperature. Store at room temperature
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10216 available at: https://www.bio-rad-antibodies.com/SDS/BUF036A
Regulatory	For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M442609:250528'

