

Datasheet: BUF037A

BATCH NUMBER 172407

Description:	ELISA GENERAL ASSAY DILUENT
Name:	ELISA GENERAL ASSAY DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Ready to use - liquid

Buffer Solution Phosphate buffered saline.

Preservative Stabilisers 0.17% ProClin 300
<10% EDTA

Intended Use

BUF037A is an ELISA assay diluent designed to reduce background noise and assist in matrix equalization in problematic serum, plasma or cell culture samples. Designed for use in all sandwich ELISA formats, it contains mammalian serum proteins to reduce non-specific interactions between sample proteins and the plate surface.

The difference in complexity, between serum or plasma samples and diluted standard curve proteins, is a potential cause of target analyte under-recovery. BUF037A equalizes the matrix, inhibits complement and thrombin activity and contains an antimicrobial agent.

Instructions For Use

Pipette buffer directly onto the plate, just prior to adding the sample or standard (suggested volume 50 - 100 ul/well).

NB. Since this buffer contains proteins derived from goat, BUF037A is not recommended

for ELISA where an anti goat secondary antibody is used.

Storage	This product is shipped at ambient temperature. Store at +4°C. DO NOT FREEZE
Guarantee	Guaranteed until date of expiry. Please see product label.
Acknowledgements	ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow.
Health And Safety Information	Material Safety Datasheet documentation #10375 available at: https://www.bio-rad-antibodies.com/SDS/BUF037A
Regulatory	For research purposes only.

Related Products

Recommended Useful Reagents

[5x ELISA COATING BUFFER \(BUF030A\)](#)

[10x ELISA WASH BUFFER \(BUF031A\)](#)

[ELISA BSA BLOCK \(BUF032A\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M439296:250523'

Printed on 29 Jan 2026