

Datasheet: BUF037C

BATCH NUMBER 158198

| | |
|----------------------|-----------------------------|
| Description: | ELISA GENERAL ASSAY DILUENT |
| Name: | ELISA GENERAL ASSAY DILUENT |
| Format: | Ready To Use |
| Product Type: | Accessory Reagent |
| Quantity: | 1000 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | Neat |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Ready to use - liquid

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.17% Proclin™ 300
<10% EDTA

Intended Use BUF037C is an ELISA assay diluent designed to reduce background noise and assist in matrix equalization in problematic serum, plasma or cell culture samples. Designed for use in all sandwich ELISA formats, it contains mammalian serum proteins to reduce non-specific interactions between sample proteins and the plate surface.

The difference in complexity, between serum or plasma samples and diluted standard curve proteins, is a potential cause of target analyte under-recovery. BUF037C equalizes the matrix, inhibits complement and thrombin activity and contains an antimicrobial agent.

Instructions For Use Pipette buffer directly onto the plate, just prior to adding the sample and standard (suggested volume 50 - 100 ul/well).

NB. Since this buffer contains proteins derived from goat, BUF037C is not recommended

for ELISA where an anti goat secondary antibody is used.

Storage

Store at +4°C.

DO NOT FREEZE

Guarantee

Guaranteed until date of expiry. Please see product label.

Acknowledgements

Proclin™ 300 is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow.

Health And Safety Information

Material Safety Datasheet documentation #10375 available at:
10375: <https://www.bio-rad-antibodies.com/uploads/MSDS/10375.pdf>

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[5x ELISA COATING BUFFER \(BUF030A\)](#)

[10x ELISA WASH BUFFER \(BUF031A\)](#)

[ELISA BSA BLOCK \(BUF032A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M379775:210412'

Printed on 08 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)