

Datasheet: BUF048

Description:	OVALBUMIN:HRP
Name:	OVALBUMIN
Format:	HRP
Product Type:	Accessory Reagent
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/50 - 1/150

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Purified Ovalbumin conjugated to Horseradish Peroxidase (HRP) - liquid
---------------------	--

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	HRP stabiliser (BUF052A) 0.01% Thiomersal
---------------------------------	--

Approx. Protein Concentrations	1.0mg/ml
---------------------------------------	----------

Intended Use	BUF048 is a preparation of HRP conjugated ovalbumin, which may be used for the measurement of ovalbumin specific IgE in immunoassay systems, without species restriction.
---------------------	---

References

- Hansen, J.S. *et al.* (2016) Limonene and its ozone-initiated reaction products attenuate allergic lung inflammation in mice. [J Immunotoxicol. 13 \(6\): 793-803.](#)
- Özkan, M. *et al.* (2021) Macrophage and dendritic cell subset composition can distinguish endotypes in adjuvant-induced asthma mouse models. [PLoS One. 16 \(6\): e0250533.](#)
- Ghonim, M.A. *et al.* (2018) Sulfated non-anticoagulant heparin blocks Th2-induced asthma by modulating the IL-4/signal transducer and activator of transcription 6/Janus

kinase 1 pathway. [J Transl Med. 16 \(1\): 243.](#)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10287 available at:
<https://www.bio-rad-antibodies.com/SDS/BUF048>
10287

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M364421:200529'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)