

Datasheet: BUF050A

BATCH NUMBER 156589

Description:	TC-PROTECTOR CELL FREEZING MEDIUM
Name:	CELL FREEZING MEDIUM
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Tissue Culture	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Cell freezing media - liquid

Preparation Chemically defined cryo-preservation medium. Does not contain FBS (foetal bovine serum) or BSA (bovine serum albumin).

Product Information **TC-Protector cell freezing medium** is a ready-to-use optimised cryo-preservation medium which allows for the direct freezing of an expansive range of mammalian cell lines at -80°C.

TC-Protector cell freezing medium provides a no-mess solution for the protection of cells during frozen storage, and reliably gives a high cell viability and durability rate on thawing.

- Instructions For Use**
1. Centrifuge cells pre-suspended in complete growth medium at 1,000-1,500rpm for 60 seconds.
 2. Remove supernatant and resuspend cells in TC-Protector at a concentration of 5×10^6 to 1×10^7 cells/mL.
 3. Transfer this cell suspension to cryogenic storage vial.

4. Freeze vial at -80°C (liquid nitrogen should be used for long-term storage).

References 1. Takebe, T. *et al.* (2014) Transient vascularization of transplanted human adult-derived progenitors promotes self-organizing cartilage. [J Clin Invest. 124 \(10\): 4325-34.](#)

Storage Store unused product at +4°C.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10150 available at: <https://www.bio-rad-antibodies.com/SDS/BUF050A>
10150

This product contains 10% dimethyl sulfoxide (DMSO): a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350304:190307'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)