

## Datasheet: BUF052A

**BATCH NUMBER 160331**

<b>Description:</b>	AbGUARD® HRP STABILIZER PLUS
<b>Name:</b>	HRP-STABILIZING DILUENT
<b>Format:</b>	Ready To Use
<b>Product Type:</b>	Accessory Reagent
<b>Quantity:</b>	100 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Product Form

Ready to use - liquid

#### Preparation

Weakly white liquid based on a Tris buffer with proprietary stabilizing components. Does not contain azide, mercury, BSA or other bovine material. Does contain weak concentrations of rabbit antibodies.

#### Preservative Stabilisers

<15ppm CMIT/MIT

#### Intended Use

AbGuard® HRP Stabilizer Plus can be used for the direct dilution and long-term stabilization/storage of Horseradish Peroxidase (HRP)-conjugated proteins and antibodies at both low and high protein concentrations.

HRP-conjugated proteins diluted with AbGuard® HRP Stabilizer Plus can be stored at 4°C for a minimum of 30 months, and can withstand temperatures up to 37°C for 2 weeks or more.

AbGuard® HRP Stabilizer Plus stabilizes the HRP enzyme and optimizes the signal-to-noise ratio in most assays. Bio-Rad recommends the use of our TMB range for this purpose.

Please note - this product cannot be used for the dilution of anti-rabbit immunoglobulin antibodies e.g Goat anti-Rabbit Ig etc.

---

<b>References</b>	1. Gholami, M. <i>et al.</i> (2015) Selenium effect on ischemia-reperfusion injury of gastrocnemius muscle in adult rats. <a href="#">Biol Trace Elem Res. 164 (2): 205-11.</a> 2. Gholami, M.R. <i>et al.</i> (2020) Protective effects of honey, <i>Apis mellifera meda.</i> skorikov, on ischemia-reperfusion induced muscle injury. <a href="#">Int J Morphol., 38 (3): 804-10,</a>
<b>Storage</b>	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10232 available at: <a href="https://www.bio-rad-antibodies.com/SDS/BUF052A">https://www.bio-rad-antibodies.com/SDS/BUF052A</a> 10232
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M360932:200127'

**Printed on 26 Jun 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)