

Datasheet: BUF053A

Description:	AbGUARD® AP STABILIZER
Name:	AP STABILIZING DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Ready to use - liquid
Preparation	Weakly white liquid based on a mildly acidic Tris buffer with proprietary components. Does not contain azide, mercury, BSA or other bovine serum material.
Preservative Stabilisers	<0.0015% CMIT/MIT
Intended Use	<p>AbGuard AP Stabilizer can be used for the direct dilution and long-term stabilization/storage of Alkaline Phosphatase (AP)-conjugated proteins and antibodies at both low and high protein concentrations.</p> <p>AP-conjugated proteins diluted with AbGuard AP Stabilizer can be stored at 4°C for a minimum of 18 months, and can withstand temperatures up to 37°C for 2 weeks or more. Avoid using phosphate buffers in assays with AP-conjugated antibodies.</p>
Storage	<p>This product is shipped at ambient temperature.</p> <p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted.</p>

Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10339 available at: https://www.bio-rad-antibodies.com/SDS/BUF053A 10339
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M439315:250523'

Printed on 23 May 2025