

Datasheet: BUF068A

Description:	AbGUARD® ANTIBODY BOOST
Name:	ANTIBODY DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Ready to use - liquid

Preservative Stabilisers None present

Intended Use AbGUARD ANTIBODY BOOST can be used to dilute concentrated or conjugated antibodies. It may be used for polyclonal or monoclonals antibodies. AbGUARD ANTIBODY BOOST amplifies the specific antibody activity up to 50-100%, with no increase in non-specific binding. For best results Bio-Rad recommends the use of our TMB substrate range. Please note - this product cannot be used for the dilution of anti-rabbit immunoglobulins e.g. Goat anti-Rabbit Ig, etc.

Instructions For Use Use a series of dilutions to determine the optimal dilution factor for individual antibodies or antibody conjugates prior to use. Once the optimal dilution factor has been determined, dilute the antibody and run the assay as normal.

Storage Store at +4°C.
DO NOT FREEZE.

Avoid direct exposure to light.

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10502 available at: https://www.bio-rad-antibodies.com/SDS/BUF068A 10502
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M407395:221007'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)