



Datasheet: BUF069A

BATCH NUMBER 155352

Description:	PROTEIN STABILIZER
Name:	PROTEIN STABILIZER
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Ready to use - liquid

Preparation Clear to slightly cloudy liquid based on a mild Tris buffer with proprietary stabilizing components. Does not contain azide, mercury or other toxic components. Does not contain bovine serum albumin or any other material derived from bovine serum.

Intended Use Protein Stabilizer stabilizes antibodies and other proteins or biomolecules in solution. It may be used to stabilize ready-to-use pre-diluted standards and controls used in immunoassays.

Instructions For Use Determine optimal dilution factors by preparing a series of dilutions of the protein in Protein Stabilizer.

Proteins diluted in Protein Stabilizer may be stored for up to 24 months at 2 - 8°C without loss of activity.

Storage Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light.

Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10519 available at: https://www.bio-rad-antibodies.com/SDS/BUF069A 10519
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350322:190307'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)