

## Datasheet: BUF069A

**BATCH NUMBER 157426**

<b>Description:</b>	PROTEIN STABILIZER
<b>Name:</b>	PROTEIN STABILIZER
<b>Format:</b>	Ready To Use
<b>Product Type:</b>	Accessory Reagent
<b>Quantity:</b>	100 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**Product Form** Ready to use - liquid

**Preparation** Clear to slightly cloudy liquid based on a mild Tris buffer with proprietary stabilizing components. Does not contain azide, mercury or other toxic components. Does not contain bovine serum albumin or any other material derived from bovine serum.

**Intended Use** Protein Stabilizer stabilizes antibodies and other proteins or biomolecules in solution. It may be used to stabilize ready-to-use pre-diluted standards and controls used in immunoassays.

**Instructions For Use** Determine optimal dilution factors by preparing a series of dilutions of the protein in Protein Stabilizer.

Proteins diluted in Protein Stabilizer may be stored for up to 24 months at 2 - 8°C without loss of activity.

**Storage** Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light.

---

**Guarantee**                      Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information**                      Material Safety Datasheet documentation #10519 available at:  
<https://www.bio-rad-antibodies.com/SDS/BUF069A>  
10519

---

**Regulatory**                      For research purposes only

---

**North & South America**    Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M350322:190307'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)