

## Datasheet: BUF069A

**BATCH NUMBER 173419**

<b>Description:</b>	PROTEIN STABILIZER
<b>Name:</b>	PROTEIN STABILIZER
<b>Format:</b>	Ready To Use
<b>Product Type:</b>	Accessory Reagent
<b>Quantity:</b>	100 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Ready to use - liquid
---------------------	-----------------------

<b>Preparation</b>	Clear to slightly cloudy liquid based on a mild Tris buffer with proprietary components. Does not contain azide, mercury or other toxic components. Does not contain bovine serum albumin or any other material derived from bovine serum.
--------------------	--

<b>Preservative Stabilisers</b>	<0.0015% CMIT/MIT
---------------------------------	-------------------

<b>Intended Use</b>	Protein Stabilizer stabilizes antibodies and other proteins or biomolecules in solution. It may be used to stabilize ready-to-use pre-diluted standards and controls used in immunoassays.
---------------------	--

<b>Instructions For Use</b>	Determine optimal dilution factors by preparing a series of dilutions of the protein in Protein Stabilizer.
-----------------------------	---

Proteins diluted in Protein Stabilizer may be stored for up to 24 months at 2 - 8°C without loss of activity.

<b>Storage</b>	This product is shipped at ambient temperature. Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10519 available at: <a href="https://www.bio-rad-antibodies.com/SDS/BUF069A">https://www.bio-rad-antibodies.com/SDS/BUF069A</a>
<b>Regulatory</b>	For research purposes only

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M450081:260304'

**Printed on 19 Mar 2026**