

Datasheet: BUF069A

Description:	PROTEIN STABILIZER
Name:	PROTEIN STABILIZER
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

Ready to use - liquid

Preparation

Clear to slightly cloudy liquid based on a mild Tris buffer with proprietary stabilizing components. Does not contain azide, mercury or other toxic components. Does not contain bovine serum albumin or any other material derived from bovine serum.

Intended Use

Protein Stabilizer stabilizes antibodies and other proteins or biomolecules in solution. It may be used to stabilize ready-to-use pre-diluted standards and controls used in immunoassays.

Instructions For Use

Determine optimal dilution factors by preparing a series of dilutions of the protein in Protein Stabilizer.

Proteins diluted in Protein Stabilizer may be stored for up to 24 months at 2 - 8°C without loss of activity.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10519 available at:
<https://www.bio-rad-antibodies.com/SDS/BUF069A>
10519

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M416145:230202'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)