

Datasheet: BUF073

BATCH NUMBER 172335

Description:	STAINING BUFFER
Name:	STAINING BUFFER
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	500 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Liquid

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Product Information **Staining buffer** has been optimized for use when staining both live and fixed cells with antibodies. It is suitable for use with both conjugated and unconjugated antibodies and can be used as an antibody dilution media and a wash buffer. The presence of bovine calf serum in this product reduces background staining to target cells. Staining buffer also contains sodium azide to minimize receptor capping.

Storage This product is shipped at ambient temperature.
Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #20370 available at: <https://www.bio-rad-antibodies.com/SDS/BUF073>

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M439333:250523'

Printed on 29 Jan 2026
