

## Datasheet: C12SAZ

<b>Description:</b>	RABBIT SERUM:Preservative Free
<b>Name:</b>	RABBIT SERUM
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Serum
<b>Quantity:</b>	20 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
Immunohistology - Resin	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Whole serum - liquid
<b>Preparation</b>	Normal rabbit serum, 0.45um filtered
<b>Preservative Stabilisers</b>	None present

**Product Information** Normal rabbit serum specifically for use in immunohistopathology. This serum has been processed to minimize haemolysis and filtered to 0.45µm to remove protein aggregates.

C12SAZ can be used as blocking sera, control sera or a serum additive in dilution buffers.

#### References

- Atalay, M.A. *et al.* (2017) Can serums be replaced by Mueller-Hinton agar in germ tube test? [Niger J Clin Pract. 20 \(1\): 61-63.](#)
- Weyrauch, A.K. *et al.* (2017) Purification of rabbit serum histidine-proline-rich glycoprotein via preparative gel electrophoresis and characterization of its glycosylation patterns. [PLoS One. 12 \(9\): e0184968.](#)

**Storage** Store at -20°C only. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the product. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10217 available at:  
10217: <https://www.bio-rad-antibodies.com/uploads/MSDS/10217.pdf>

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M364454:200529'

**Printed on 29 Aug 2021**

---

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)