

Datasheet: C13SD

BATCH NUMBER 165472

Description:	RAT SERUM WITH AZIDE
Name:	RAT SERUM
Format:	Serum
Product Type:	Serum
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Whole serum - liquid
---------------------	----------------------

Preparation	Normal rat serum.
--------------------	-------------------

Preservative Stabilisers	0.1% Sodium Azide
---------------------------------	-------------------

References

1. Busche, A. *et al.* (2011) Genetic labeling reveals altered turnover and stability of innate lymphocytes in latent mouse cytomegalovirus infection. [J Immunol. 186: 2918-25.](#)
2. Dorrell, C. *et al.* (2011) Prospective isolation of a bipotential clonogenic liver progenitor cell in adult mice. [Genes Dev. 25: 1193-203.](#)
3. Chennupati, V. *et al.* (2010) Intra- and intercompartmental movement of gammadelta T cells: intestinal intraepithelial and peripheral gammadelta T cells represent exclusive nonoverlapping populations with distinct migration characteristics. [J Immunol. 185: 5160-8.](#)
4. Gelosa, P. *et al.* (2010) Peroxisome proliferator-activated receptor {alpha} agonism prevents renal damage and the oxidative stress and inflammatory processes affecting the brains of stroke-prone rats. [J Pharmacol Exp Ther. 335: 324-31.](#)

5. Wilson, S. *et al.* (2009) Hepatosplenomegaly associated with chronic malaria exposure: evidence for a pro-inflammatory mechanism exacerbated by schistosomiasis. [Parasite Immunol. 31: 64-71.](#)
6. Muschaweck, M. *et al.* (2021) Cognate recognition of microbial antigens defines constricted CD4⁺ T cell receptor repertoires in the inflamed colon. [Immunity. S1074-7613\(21\)00346-0.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the product. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10496 available at: <https://www.bio-rad-antibodies.com/SDS/C13SD10496>

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M417465:230330'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)