

Datasheet: FCSC815B

BATCH NUMBER 167036

Description:	QUANTUM™ SIMPLY CELLULAR® MOUSE IgG
Name:	QUANTUM™ SIMPLY CELLULAR® IgG
Format:	Flow Cytometry Calibration Reagent
Product Type:	Accessory Reagent
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Buffer Solution	Borate buffered saline.
------------------------	-------------------------

Preservative	0.1% Bovine Serum Albumin
Stabilisers	0.05% Tween 20
	10mM EDTA
	0.09% Sodium Azide (NaN ₃)

Product Information **Quantum™ Simply Cellular® Mouse IgG** is comprised of four coated populations and one blank population of uniform microspheres that are approximately the size of human lymphocytes (7-9µm). Each of the four sets is coated with goat-anti-mouse IgG (Fc-specific) with different calibrated Antibody Binding Capacities (ABC) for mouse monoclonal antibodies. The blank population has no specific binding capacity for mouse IgG.

When the bead populations are labeled in the same manner as the cells to be analyzed, they provide a means for constructing a QuickCal® calibration curve (ABC values vs. fluorescence intensity), from which samples may be "read". Cellular antigen expression may thus be quantitated in ABC units. The calibration curve is also a useful means for determining instrument linearity.

Intended Use FCSC815B is suitable for the quantitation of cellular antigen expression in Antibody Binding Capacity (ABC) units. This product may also serve as an accurate compensation standard when labeled with the same antibodies used to stain cell samples.

Reagents In The Kit 1 x 5ml bottle of uncoated microbeads
4 x 5ml bottles of coated microbeads.

Instructions For Use Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/FCSC815B.pdf

References

1. Park, J.A. & Cheung, N.V. (2022) Overcoming tumor heterogeneity by *ex vivo* arming of T cells using multiple bispecific antibodies. [J Immunother Cancer. 10 \(1\) :e003771.](https://doi.org/10.1186/s12943-022-00377-1)
2. Bryniarski, M.A. *et al.* (2024) Utility of Cellular Measurements of Non-Specific Endocytosis to Assess the Target-Independent Clearance of Monoclonal Antibodies [bioRxiv 16 Apr \[Preprint\].](https://doi.org/10.1101/2024.04.16.596111)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. This product is photosensitive and should be protected from light.

Guarantee Guaranteed until date of expiry. Please see product label.

Acknowledgements Quantum and Simply Cellular are trademarks of Bangs Laboratories, INC.

Health And Safety Information Material Safety Datasheet documentation #10256 available at: <https://www.bio-rad-antibodies.com/SDS/FCSC815B>
10256

Regulatory For research purposes only.

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M416244:230220'

Printed on 16 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)