

Datasheet: HCA037

BATCH NUMBER 150051

Description:	HUMAN ANTI MOUSE IgG2a
Specificity:	IgG2a
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD06344
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/5000 - 1/12500
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	A bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library. Expressed in <i>E. coli</i> and purified using NiNTA affinity chromatography. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal

Approx. Protein Concentrations	0.5 mg/ml
Immunogen	Purified mouse monoclonal IgG2a.
External Database Links	<p>UniProt:</p> <p>P01865 Related reagents</p> <p>P01863 Related reagents</p> <p>P01864 Related reagents</p> <p>Entrez Gene:</p> <p>380793 Igh-1a Related reagents</p> <p>380793 Igh-1a Related reagents</p> <p>380793 Igh-1a Related reagents</p>
RRID	AB_906057
Specificity	Human anti Mouse IgG2a, clone AbD06344 is specific for mouse IgG2a. No cross reactivity is seen with other mouse immunoglobulin sub-classes, rat immunoglobulins, or sheep, bovine and goat IgG.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.
Guarantee	6 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA03710094
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Mouse Anti Human C-MYC (MCA2200...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M371478:200612'

Printed on 06 Apr 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)