

Datasheet: HCA050A

Description:	RECOMBINANT HUMAN IgG4 LAMBDA
Name:	HUMAN IgG4 LAMBDA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD00264_hIgG4
Isotype:	IgG4 Lambda
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG4 antibody selected from the HuCAL GOLD phage display library and expressed in a human cell line, HKB-11 - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml

Immunogen	Green fluorescent protein
External Database Links	<p>UniProt: P01861 Related reagents</p> <p>Entrez Gene: 3503 IGHG4 Related reagents</p>
RRID	AB_915465
Product Information	<p>Recombinant Human IgG4 lambda, clone AbD00264_hlgG4 is a recombinant human IgG4 antibody with a lambda light chain.</p> <p>Clone AbD00264_hlgG4 is specific for green fluorescent protein (GFP) and has no known reactivity with mammalian proteins or other antigens. This product is recommended for use as a standard in assays designed to measure IgG4 levels, or as a control antibody when using other human antibodies of the same isotype and subclass.</p> <p>In addition to the tetrameric IgG structure, composed of two copies each of the heavy and light chains, human IgG4 forms so-called half-molecules. This species contains only one heavy and one light chain. Half-molecules can assemble with an unrelated half-molecule to form bispecific tetrameric antibody species (Aalberse et al. 2002). Half-molecules of human IgG4 preparations can be detected as an additional species under denaturing, non-reducing conditions using SDS-PAGE, size exclusion chromatography or capillary electrophoresis.</p> <p>Mutation in the core hinge region of IgG4 prevents the formation of half-molecules (Bloom et al. 1997). In addition to the wild type recombinant, Bio-Rad also offers recombinant Human IgG4 lambda in its mutated form (HCA246). Recombinant Human IgG4 kappa is also available in both wild type (HCA195) and mutant form (HCA247).</p>
Further Reading	<ol style="list-style-type: none"> 1. Aalberse, R.C. & Schuurman, J. (2002) IgG4 breaking the rules. Immunology. 105 (1): 9-19. 2. Bloom, J.W. et al. (1997) Intrachain disulfide bond in the core hinge region of human IgG4. Protein Sci. 6:407-415
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See

bio-rad.com/en-us/trademarks for details.

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA050A 10040
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Secondary Antibodies

Goat Anti Human IgG F(ab')₂ (0500-0099...) [HRP](#)

Mouse Anti Human IgG (Fc) (CH2 DOMAIN) (MCA647...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M428274:240301'

Printed on 07 Apr 2025