

Datasheet: HCA160

BATCH NUMBER 172175

Description:	HUMAN ANTI PLASMODIUM FALCIPARUM HRPII
Specificity:	HISTIDINE RICH PROTEIN II
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD14971_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Protozoan
Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This product is supplied as liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% sodium azide (NaN ₃)

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant HRPII fusion protein from <i>Plasmodium falciparum</i> .
External Database Links	UniProt: Q25878 Related reagents
RRID	AB_10896836
Specificity	Human anti <i>Plasmodium falciparum</i> HRPII, clone AbD14971 recognizes HRPII (histidine rich protein II) from the Malaria causing protozoan species <i>Plasmodium falciparum</i> . Clone AbD14971 was selected by phage display under heat stress conditions. The T _m value for thermal denaturation of the antibody was measured by CD spectroscopy and has been found to be 70 °C.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA160
Licensed Use	For <i>in vitro</i> . research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M428303:240301'

Printed on 06 Mar 2026

