

Datasheet: HCA166

BATCH NUMBER 148925B

Description:	HUMAN ANTI TRASTUZUMAB
Specificity:	TRASTUZUMAB
Other names:	HERCEPTIN®
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD16712
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library, expressed in E. coli. The antibody is tagged with a V5 tag and a double extended Strep-tag at the C-terminus of the antibody heavy chain. Supplied as a liquid.

Preparation

StrepTactin affinity chromatography

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Herceptin

Specificity	<p>Human anti Trastuzumab antibody, clone AbD16712 a monoclonal antibody recognizes Trastuzumab, also known as Herceptin, is a drug used in the treatment of HER2 positive breast cancer and other HER2 overexpressing cancers including HER2-positive metastatic cancers of the gastrointestinal tract. Herceptin binds to the HER2 (or c-erbB2) proto-oncogene, an EGF receptor-like protein found on 20-30% of breast cancer cells. The binding leads to antibody mediated (complement mediated) killing of the HER2 positive cells.</p> <p>Human anti Trastuzumab (AbD16712) is an anti-idiotypic antibody that specifically recognizes the Herceptin humanised monoclonal antibody and can be used to measure the Heceptin levels in serum from patients.</p> <p>This clone has been demonstrated to inhibit the Herceptin - HER2 interaction. Removal of the preservative is recommended prior to use in functional assays.</p>
Affinity	The monovalent intrinsic affinity of this antibody was measured as $KD=0.2$ nM by real time, label-free molecular interaction analysis on immobilized trastuzumab.
ELISA	This product may be used as a detection antibody, or as capture antibody in a sandwich ELISA together with HCA168 , when conjugated to HRP, as the detection reagent. Protocol: PK bridging ELISA to measure free drug
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Herceptin® is a registered trademark of Genentech, Inc.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA166 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG IMMO (MCA2488...) [Purified](#)

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Biotin](#), [HRP](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[HUMAN ANTI TRASTUZUMAB \(HCA167\)](#)

[HUMAN ANTI TRASTUZUMAB \(HCA168\)](#)

[HUMAN ANTI TRASTUZUMAB \(HCA169\)](#)

[HUMAN ANTI TRASTUZUMAB \(HCA176\)](#)

[HUMAN ANTI TRASTUZUMAB \(HCA177\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M372729:200715'

Printed on 01 Mar 2024
