

Datasheet: HCA171G

BATCH NUMBER 163112

Description:	RECOMBINANT HUMAN IgE LAMBDA
Name:	HUMAN IgE LAMBDA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD00264_hIgE
Isotype:	IgE Lambda
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgE antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgE prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 1.0 mg/ml

Immunogen Green fluorescent protein

Product Information **Recombinant Human IgE lambda, clone AbD00264_hlgE** contains a lambda light chain.

Clone AbD00264_hlgE is a recombinant antibody with specificity for green fluorescent protein (GFP). It has no known reactivity with mammalian proteins or other antigens. It is therefore recommended as a control reagent for assays where this immunoglobulin is measured, or when using other human antibodies of the same isotype and subclass. It is recommended that this reagent is used at the same concentration as the test reagent.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Health And Safety Information Material Safety Datasheet documentation #10094 available at: <https://www.bio-rad-antibodies.com/SDS/HCA171G10094>

Licensed Use For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391646:211015'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)