

Datasheet: HCA171G

BATCH NUMBER 174467

Description:	RECOMBINANT HUMAN IgE LAMBDA
Name:	HUMAN IgE LAMBDA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD00264_hIgE
Isotype:	IgE Lambda
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgE antibody selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgE prepared by affinity chromatography
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein	Antibody concentration 1.0 mg/ml

Concentrations

Immunogen	Green fluorescent protein
------------------	---------------------------

Product Information	Recombinant Human IgE lambda, clone AbD00264_hIgE contains a lambda light chain.
----------------------------	---

Clone AbD00264_hIgE is a recombinant antibody with specificity for green fluorescent protein (GFP). It has no known reactivity with mammalian proteins or other antigens. It is therefore recommended as a control reagent for assays where this immunoglobulin is measured, or when using other human antibodies of the same isotype and subclass. It is recommended that this reagent is used at the same concentration as the test reagent.

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
-------------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA171G
--------------------------------------	--

Licensed Use	For <i>in vitro</i> . research purposes only, unless otherwise specified in writing by Bio-Rad.
---------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .
-------------------------	---

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M450361:260312'

Printed on 24 Aug 2026