

Datasheet: HCA173F

BATCH NUMBER 1709

Description:	HUMAN ANTI CHICKEN CD25:FITC
Specificity:	CD25
Other names:	IL-2 RECEPTOR ALPHA
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	AbD13504
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10
Immunofluorescence			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Chicken

Species Cross Reactivity Does not react with:Duck, Goose

Product Form A bivalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in *E. coli*. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. Conjugated to fluorescein isothiocyanate (FITC) - liquid.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

Preparation Metal chelate affinity chromatography

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	Antibody concentration 0.1 mg/ml
Immunogen	Recombinantly expressed his-tagged chicken CD25
External Database Links	UniProt: Q90WJ4 Related reagents
Specificity	Human anti Chicken CD25 antibody, clone AbD13504 recognizes the chicken IL-2 receptor alpha chain (IL-2R α) , identified as the chicken homologue of human CD25 and cited as ChCD25 in published literature. Chicken IL-2R α is a cell surface antigen expressed by activated T-lymphocytes, monocytes/macrophages and thrombocytes. CD4+ve/CD25+ve and CD8+ve/CD25+ve cells are up-regulated following infection with avian influenza virus (H9N2) and may therefore have an important role in the immune response to avian influenza virus infection (Huang et al. 2012). Chicken IL-2R α is the low affinity receptor for Interleukin-2.
References	1. Braukmann, M. <i>et al.</i> (2015) Avian CD25(+) gamma/delta ($\gamma\delta$) T cells after <i>Salmonella</i> exposure. Vet Immunol Immunopathol. 168 (1-2): 14-8. 2. Polasky C <i>et al.</i> (2016) Non-specific activation of CD8 α -characterised $\gamma\delta$ T cells in PBL cultures of different chicken lines. Vet Immunol Immunopathol. 179: 1-7.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/HCA173F10041
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[MOUSE ANTI CHICKEN CD4:RPE \(MCA2164PE\)](#)

[MOUSE ANTI CHICKEN CD8 ALPHA:RPE \(MCA2166PE\)](#)

[MOUSE ANTI CHICKEN CD25 \(MCA5925GA\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371577:200612'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)