

Datasheet: HCA181

Description:	HUMAN ANTI HUMAN HbA1c
Specificity:	HbA1c
Other names:	Glycated hemoglobin, Hb1c, Hemoglobin A1c
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD15790
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in *E. coli*. The antibody is tagged with a V5 tag and a double extended Strep-tag at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

StrepTactin affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Total protein concentration 0.5 mg/ml

Immunogen	The immunogens used in the generation of this antibody were: A glycosylated peptide, corresponding to N-terminal 9 amino acids of haemoglobin beta chain coupled to carrier proteins Purified HbA1c
External Database Links	UniProt: P68871 Related reagents Entrez Gene: 3043 HBB Related reagents
Specificity	Human anti Human HbA1c antibody, clone AbD15790 recognizes HbA1c, i.e. hemoglobin when glycosylated at the N-terminal valine of the hemoglobin β chain. It does not bind to non-glycosylated hemoglobin. Besides binding to the normal glycosylated hemoglobin, it also recognizes the glycosylated hemoglobin variants HbE and HbD, but does not recognize the HbS or HbC variants. Elevated HbF cannot be detected due to the absence of the hemoglobin β chain. The measurement of HbA1c levels is used to identify the average plasma glucose concentration. In Diabetes mellitus, higher levels of HbA1c have been associated with cardiovascular disease and other diseases. Levels of more than 6.5% HbA1c (48 mmol/mol) are used as a criterion to diagnose diabetes.
Affinity	The affinity of this antibody, as measured against the immunogen, was measured as $K_D=13$ nM by real time, label-free molecular interaction analysis using HbA1c-BSA and HbA1c-TRF conjugates as immobilized antigen. The affinity of this antibody, as measured against purified HbA1c, was measured as $K_D=80$ nM by real time, label-free molecular interaction analysis using full length HbA1c as immobilized antigen.
ELISA	This product is suitable for use indirect or indirect ELISA applications
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at $2-8^\circ\text{C}$ for short term use (up to 4 weeks) and store the remaining aliquots at -20°C . Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety	Material Safety Datasheet documentation #10094 available at:

Information	10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Alk. Phos.](#), [Biotin](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391663:211015'

Printed on 15 Oct 2021